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
The screenshot shows a web-based interface for managing DICOM data. At the top, there is a toolbar with various icons. Below it, the title 'Base de datos DICOM' is displayed. The main area contains a search bar with fields for 'Pacientes', 'Estudios', and 'Series'. Below the search bar is a table with the following data:

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Below the table, there is a large orange text box with the following text: 'Dar clic en Aceptar para subir el directorio en la base de datos DICOM'. At the bottom right of the interface, the 'Aceptar' button is circled in red. Below the table, there is a status bar that reads 'Importación efectuada: añadidos 1 pacientes, 1 estudios, 1 series, 291 instancias.' Below this, there are tabs for 'Datos DICOM' and 'Lector', and a section for 'Advertencias'. At the bottom of the interface, there are buttons for 'Desmarcar todo', 'Examinar', and 'Cargar', along with a checkbox for 'Avanzado' which is checked. A copyright notice '© Sonia Pujol, Ph.D. All Rights Reserved' is visible at the bottom center, and a page number '10' is at the bottom right.

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The screenshot displays the 3D Slicer 5.6.2 software interface. The top menu bar includes 'Archivo', 'Editar', 'Vista', and 'Ayuda'. Below the menu is a toolbar with various icons for file operations and viewing. The main workspace is divided into several panels:


- Left Panel (Data and Settings):**
 - 'Ayuda y reconocimiento' section.
 - 'Datos cargados' (Loaded Data) section showing a tree view with 'patient1 (patient1_ID)', 'CT Thorax Abdomen (20050601)', and '6: CT_Thorax_Abdomen'.
 - 'Conexión a red DICOM' (DICOM Network Connection) section.
 - 'Configuración de la base de datos DICOM' (DICOM Database Configuration) section.
 - 'Plugins DICOM' section.
 - 'Coordenadas' (Coordinates) section with 'Mostrar corte ampliado' (Show enlarged slice) checkbox and 'L', 'F', 'B' labels.
- Top Right Panel (3D View):** A 3D perspective view of a purple rectangular volume with axes labeled 'S' (Superior), 'I' (Inferior), 'R' (Right), and 'L' (Left). A 'P' (Posterior) label is also visible.
- Bottom Panel (MPR):** Three orthogonal views of the CT scan:
 - Left (Axial):** Labeled 'R' (Right) and 'S' (Superior). Position: 'SI -188.500mm'.
 - Middle (Coronal):** Labeled 'G' (Grayscale) and 'A' (Anterior). Position: 'A: 169.2539mm'.
 - Right (Sagittal):** Labeled 'Y' (Yellow) and 'L' (Left). Position: 'L: -6.7461mm'.

A green text box on the left side of the interface contains the following text:

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The image shows the 3D Slicer 5.6.2 software interface. The main window displays a 3D volume rendering of a CT scan, with a pink wireframe box indicating the current slice position. The volume is labeled with 'S' for Superior, 'R' for Right, and 'L' for Left. Below the 3D view are three 2D slice views: Axial (left), Coronal (middle), and Sagittal (right). Each slice view has its own window with a title bar and a zoom slider. The Axial view is titled '6: CT_Thorax_Abdomen' and shows a cross-section of the abdomen. The Coronal view is titled '6: CT_T...Abdomen' and shows a front view of the abdomen. The Sagittal view is titled '6: CT_T...Abdomen' and shows a side view of the abdomen. The interface includes a menu bar (Archivo, Editar, Vista, Ayuda), a toolbar with various icons, and a left sidebar with several panels: 'Ayuda y reconocimiento', 'Visualización' (with a filter and search table), 'Histograma', 'Leyenda de colores', and 'Coordenadas'. A green text box is overlaid on the left sidebar, containing the text: 'Los tres segmentos anatómicos aparecen en el visor 3D.'

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


The screenshot displays the 3D Slicer 5.6.2 interface. The main 3D viewer shows a volume rendering of a CT scan with a purple bounding box and orientation labels: 'S' (Superior) at the top, 'I' (Inferior) at the bottom, 'R' (Right) on the left, and 'L' (Left) on the right. To the right, a multi-planar reconstruction (MPR) view shows three orthogonal slices: Axial (top), Coronal (middle), and Sagittal (bottom). The software interface includes a menu bar (Archivo, Editar, Vista, Ayuda), a toolbar with various icons, and a left sidebar with panels for 'Filtros', 'Tabla de búsquedas', 'Interpolación', 'Ventana/nivel', 'Umbral', 'Histograma', 'Leyenda de colores', and 'Coordenadas'. The 'Ventana/nivel' panel shows a manual window level of 350 and a width of 40. The 'Umbral' panel shows a threshold of -1024 to 3071.

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

The image shows the 3D Slicer 5.6.2 interface. On the left, the 'Renderización de volumen' (Volume Rendering) tab is active. The 'Filtro' (Filter) section shows a grid of rendering modes. The 'CT-Cardiac3' filter is selected, highlighted with a red arrow. The 'Cortar' (Cut) section shows 'CT-Cardiac3' selected. The 'Renderiza' (Render) section shows 'CT-Cardiac3' selected. The 'Pará' (Parameters) section shows '8' and '16' selected. The 'Coordenadas' (Coordinates) section shows 'Mostrar corte ampliado' (Show enlarged cut) checked. The 'L', 'F', and 'B' checkboxes are also visible. In the center, a 3D view of a CT scan is shown, with a purple wireframe box around it. The axes are labeled 'A', 'P', and 'I'. On the right, a multi-panel view shows three axial CT slices of the thorax and abdomen, with a red, green, and yellow progress bar above each panel. The top panel is labeled 'Si -170.500mm', the middle 'A: 167.1354mm', and the bottom 'L: -64.1289mm'. The text 'B: 6: CT_Thorax_Abdomen' is visible in each panel.

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