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

# Carregando um volume DICOM

The screenshot shows the 3D Slicer 5.6.2 interface. The main window displays the DICOM database panel with the following data:

Pacientes:	Estudos:	Séries:			
Nome do paciente	ID do pacie	Data de nascimento	Sexo	Estudos	Data do últi
patient1	pati...1_ID			1	2005-06-01 2024...782
Gamotis, Robert	10841	1981-03-19	M	0	2023...683

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
# Carregando um volume DICOM

The screenshot displays the 3D Slicer 5.6.2 software interface. The top menu bar includes 'Arquivo', 'Editar', 'Visualizar', and 'Ajuda'. The main toolbar contains various icons for file operations and viewing. The left sidebar shows the 'Dados Carregados' (Loaded Data) section with a tree view containing 'patient1 (patient1\_ID)', 'CT Thorax Abdomen (20050601)', and '6: CT\_Thorax\_Abdomen'. Below this, there are sections for 'Rede DICOM', 'Configurações do banco de dados DICOM', 'Plugins DICOM', and 'Sonda de Dados'. The main 3D view area shows a purple background with a white wireframe box representing the volume, labeled with 'S' (Superior), 'I' (Inferior), 'R' (Right), and 'L' (Left). Below the 3D view, there are three multi-planar views: an axial view (labeled 'S: -188.5000mm'), a coronal view (labeled 'A: 169.2539mm'), and a sagittal view (labeled 'L: -6.7461mm'). Each view shows a grayscale CT scan of the thorax and abdomen. The bottom left corner of the interface has a small 'Abdomen' label.

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Selecionar o módulo **Volumes**

# Carregando um volume DICOM


The screenshot displays the 3D Slicer 5.6.2 interface. The top menu bar includes 'Arquivo', 'Editar', 'Visualizar', and 'Ajuda'. Below it is a toolbar with various icons. The left sidebar contains the following sections:

- 3D Slicer**
- Ajuda e Agradecimentos**
- Volume Ativo: 6: CT\_Thorax\_Abdomen**
- Informações do volume**
- Mostrar**
- Predefinições:** A row of six icons representing different CT scan presets. The fourth icon, showing a cross-section of the abdomen, is circled in red.
- Tabela de cores:** A slider set to 'Grey'.
- Interpolação:** A dropdown menu set to 'v'.
- Janelamento/Nível:** A slider set to 'Auto' with 'W: 2217' and 'L: 85'.
- Limiar:** A dropdown menu set to 'Desligado' with a slider ranging from '-1024' to '3071'.
- Histograma**
- Legenda de cores**
- Sonda de Dados**
- Mostrar Fatia Ampliada**
- L**
- F**
- B**

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At the bottom, three orthogonal views are displayed: a coronal view (left), a sagittal view (middle), and a transverse view (right). Each view has a corresponding color-coded slider above it: 'R' (Red) for the coronal view, 'S' (Green) for the sagittal view, and 'Y' (Yellow) for the transverse view. The sagittal view shows a value of 'A: 169.2539mm' and the transverse view shows 'L: -6.7461mm'. The labels 'B: 6: CT\_Thorax\_Abdomen' are visible at the bottom of each view.

# Carregando um volume DICOM

The screenshot displays the 3D Slicer 5.6.2 interface. The left sidebar contains the 'Volumes' panel with the active volume '6: CT\_Thorax\_Abdomen'. Below it are 'Predefinições' (Presets) and 'Tabela de cores' (Color Table) set to 'Grey'. The 'Janelamento/Nível' (Window/Level) section is set to 'Manual' with 'W: 350' and 'L: 40'. The 'Limiar' (Threshold) is set to 'Desligado' (Off) with a slider from -1024 to 3071. The main 3D view shows a CT scan slice with a red box highlighting the slice menu controls. The slice menu includes 'R' (Right), 'A' (Anterior), and 'I' (Inferior) icons, a slider for 'S: -188.5000mm', and a 'G' (Grayscale) icon. The 'Axial' view is selected. Three red arrows point to the 'R', 'A', and 'I' icons. The bottom of the interface shows three small thumbnail views of the CT scan in different orientations.

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B: 6: CT\_Thorax\_Abdomen

B: 6: CT\_Thorax\_Abdomen

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The image shows the 3D Slicer 5.6.2 software interface. On the left, there is a sidebar with various settings and a 'Módulos' (Modules) section. A large orange text box is overlaid on the left side, containing the instruction: 'Clique no ícone do menu do layout do Slicer e selecione Widescreen convencional'. The main window displays a CT scan of a thorax and abdomen. The top menu bar includes 'Arquivo', 'Editar', 'Visualizar', and 'Ajuda'. The 'Módulos' dropdown menu is open, showing a list of layout options. The 'Widescreen Convencional' option is highlighted. Below the menu, there are two smaller preview windows showing different views of the CT scan. The bottom status bar shows coordinates: 'G: 169.2539mm', 'Y: -6.7461mm', and 'L: -6.7461mm'.

3D Slicer 5.6.2

Arquivo Editar Visualizar Ajuda

Módulos: Volumes

3D Slicer

Ajuda e Ag

Volume Ativo:

Informaçõ

Mostrar

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Tabela de cores: Grey

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W: 350 L: 40

Limiar: Desligado

-1024 3071

Histograma

Legenda de cores

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Mostrar Fatia Ampliada

L

F

B

Conventional

Widescreen Convencional

Gráfico Convencional

Quatro Quadros

Quatro Quadros c/ Tabela

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Quatro-Quadros Quantitativo

3D duplo

Triplo 3D

Apenas 3D

Tabela 3D

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Apenas a fatia amarela

Apenas a fatia verde

3D com guias

Fatia com guias

Comparar

Comparar Widescreen

Comparar Grade

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Dois sobre dois

Lado a lado

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Fatias quatro por dois

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Quatro-em-um, monitor duplo


G A: 169.2539mm Y L: -6.7461mm

B: 6: CT\_Thorax\_Abdomen


B: 6: CT\_Thorax\_Abdomen

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

# Controlador do Visor 3D

The image shows the 3D Slicer software interface. The top menu bar includes 'Arquivo', 'Editar', 'Visualizar', and 'Ajuda'. Below the menu is a toolbar with various icons. The main window is divided into several panels. On the left, there is a sidebar with '3D Slicer' and 'Ajuda e Agradecimentos'. Below that, it shows 'Volume Ativo: 6: CT\_Thorax\_Abdomen' and 'Informações do volume'. The central 3D viewer shows a 3D reconstruction of a CT scan of the thorax and abdomen. The right side of the interface features three orthogonal views: a coronal view at the top, an axial view in the middle, and a sagittal view at the bottom. Each view has a corresponding control bar with a slider and a label (S, A, L). The 3D viewer control panel is located at the top of the 3D viewer area, featuring a coordinate system (P, S, L, R, A, I) and several icons. Two red arrows point to specific icons in this panel: one points to the 'Pin' icon (a vertical line with a dot) and the other points to the 'Reset View' icon (a circular arrow).

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



The image shows the 3D Slicer software interface. The main window displays a 3D volume rendering of a CT scan, showing the thorax and abdomen. The volume is rendered in a semi-transparent brown color. A red arrow points to a slider control in the left sidebar, which is used to adjust the transfer function and display the aorta. The sidebar shows the following settings:

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The right sidebar shows the multi-planar reconstructions (MPR) of the CT scan:

- Axial: 6: CT\_Thorax\_Abdomen (S: -192.5190mm)
- Coronal: 6: CT\_Thorax\_Abdomen (A: 157.5740mm)
- Sagittal: 6: CT\_Thorax\_Abdomen (L: 5.7703mm)

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





The image shows the 3D Slicer 5.6.2 interface. On the left, the 'Módulos' (Modules) list is visible, with 'Modelos' (Models) highlighted in blue. A red arrow points to 'Modelos' in the list. The main window displays a 3D surface model of a human head in a neutral view. To the right, three axial MRI slices are shown in grayscale. The top slice is labeled 'S: -27.1875mm', the middle slice 'A: -34.5000mm', and the bottom slice 'L: -14.0625mm'. The interface includes a menu bar (Arquivo, Editar, Visualizar), a toolbar, and various control panels for visibility and display settings.

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